



## Analysis of Organotin Compounds via QuEChERS and LC-MS/MS

### - Brief Description -

#### Sample extraction (QuEChERS):

Detailed guidelines for the sample preparation (QuEChERS method) are given in the CEN method EN 15662 (Determination of pesticide residues using GC-MS and/or LC-MS (/MS) following acetonitrile extraction/partitioning and cleanup by dispersive SPE - QuEChERS method) or as brief description on the Internet ([www.quechers.de](http://www.quechers.de)).

#### Analysis by LC-MS/MS:

(Please consider: the LC-MS/MS data is exemplary)

<b>Analyte scope</b>	Cyhexatin (sum of azocylotin and cyhexatin) Fentin Fenbutatin oxide Internal Standard: TPP																		
<b>Instrument parameters</b>	API 4000																		
<b>Ionisation mode</b>	ESI pos																		
<b>Column</b>	Zorbax 3,5 µm; Eclipse XDB-C18; 2,1x 50 mm																		
<b>Pre-column</b>	C18 ODS 4mm x 2mm ID (Phenomenex AJO-4286)																		
<b>Eluent A</b>	5 mmol NH <sub>4</sub> formiate in H <sub>2</sub> O +1% formic acid																		
<b>Eluent B</b>	5 mmol NH <sub>4</sub> formiate in Methanol +1% formic acid																		
<b>Gradient</b>	<p>flow: 0.4ml/min</p> <table border="1"> <thead> <tr> <th>time</th> <th>A</th> <th>B</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>60% to</td> <td>40%</td> </tr> <tr> <td>2</td> <td>0% to</td> <td>100%</td> </tr> <tr> <td>8</td> <td>0% to</td> <td>100%</td> </tr> <tr> <td>8.1</td> <td>60% to</td> <td>40%</td> </tr> <tr> <td>13</td> <td>60% to</td> <td>40%</td> </tr> </tbody> </table>	time	A	B	0	60% to	40%	2	0% to	100%	8	0% to	100%	8.1	60% to	40%	13	60% to	40%
time	A	B																	
0	60% to	40%																	
2	0% to	100%																	
8	0% to	100%																	
8.1	60% to	40%																	
13	60% to	40%																	

## Acquired mass transitions

Analyte, mass transition (m/z)	Relative Sensitivity 6= best)	First Mass	Second Mass		DP	CE	CXP
Cyhexatin 365/201	5	365.1	200.9		71	27	30
Cyhexatin 365/283	6	365.1	283.1		71	17	8
Cyhexatin 367/203	3	367.2	202.9		71	25	28
Cyhexatin 367/285	4	367.2	285.0		71	17	10
Cyhexatin 369/205	1	369.2	205.0		71	23	18
Cyhexatin 369/287	2	369.2	287.0		71	15	18
Fentin 347/116	6	347.0	115.9		96	43	20
Fentin 347/193	5	347.0	193.0		96	37	16
Fentin 349/118	4	349.0	117.9		91	41	20
Fentin 349/195	3	349.0	195.0		91	37	16
Fentin 351/120	2	351.0	119.9		96	41	8
Fentin 351/197	1	351.0	197.1		96	37	14
Fenbutatin oxide 517/195	6	517.3	195.0		116	79	18
Fenbutatin oxide 517/349	5	517.3	349.0		116	45	12
Fenbutatin oxide 517/461	4	517.3	461.1		116	35	14
Fenbutatin oxide 519/197	2	519.3	197.0		121	67	16
Fenbutatin oxide 519/351	3	519.3	350.9		121	47	20
Fenbutatin oxide 519/463	1	519.3	463.3		121	35	14

## Quantification:

The analyte concentration is calculated via regression curve of the area ratios of the analyte to the internal standard. In case of suspected MRL-violations a standard addition procedure is used.

## Performance:

Reporting levels using API 4000 instrumentation: 0.002 mg/kg for all compounds

Validation data CVUA Stuttgart brief overview:

Pesticide	Level min	Level max	Rec Median	Rec Mean	CV [%]	# of rec	# of outl	% Rec (70-120%)	# of Labs
Cyhexatin	0,1	0,1	<b>83</b>	84	11,9	20	0	<b>90</b>	1
Fenbutatin Oxide	0,04	0,1	<b>93</b>	91	13,3	21	0	<b>95</b>	1
Fentin	0,1	0,1	<b>86</b>	86	9	20	0	<b>100</b>	1

For more data see [www.crl-pesticides-datapool.eu](http://www.crl-pesticides-datapool.eu), method validation data